

AHF-Procurement -2023-13 - ملحق دفتر الشروط / Tender Book annex				C.4-PR236 (ERCP)		
/ Commercial Name / Ticari Ad الإسم التجاري		Return Quotation to /T				
no رق م	اسم المواد ووصفها ITEMS DESCRIPTION/SPECIFICAT ION Malzemenin Adı ve Özelliği	المواصفات الفنية Specification Teknik Özellikler		الملاحظات Notes Notlar		العلامة التجارية & بلد المنشأ Brand Details & Country of Originrı
1	ERCP	- Second hand (Used device): Olympus - TJF-Q180V ERCP with operating table		جهاز مستعمل - نوع اولومبوس حسب الموديل المذكور بحيث يتوافق مع برج التنظير المتوفر في المشفى		
3	Cuttry device	Product Technical specification :Six Output working modes a) Cut: 400W b) Blend 1: 250W c) Blend 2: 200W d) Spray coagulation: 100W e) Soft coagulation: 80W f) Bipolar coagulation: 50W		الماني او تركي المنشأ ضمن المواصفات المذكورة - يستخدم لجهاز ال ERCP مزود بمدخل يتوافق مع جهاز ال ERCP + الوصلة الخاصة مابين الكوتري وجهاز ال ERCP		
		Used X-ray generat High frequency inverter frequency from the property of the voltage: Radiography tube voltage: Radiography tube voltage: Rotating anode, dual focus, small 0. Collimator Collimator Collimator Collimator Flectric lung field iris copper sheet, iris .counterclockwise/clockwise, can be Advanced technology: Liquid metal I Dynamic flat panel det Type: Amorphous silice Nominal area:210mm*2 Pixel pitch: 205µm /Matrix: C-arm gantry /Lateral movem Vertical lifting movement Swing left & right: 25°(±12°, (90°~45°-)C-arm orbital radiography	sency:40kHz 50mA、 0.1s) 40-120KV 40-120KV Bmm, large 0.6mm ID continuous adjustable collimator can rotate opened and closed bearing technology ector on CSI 10mm 1024×1024 ent:250mm error ±1°)			



Rotate around the horizontal axis: ±225° LAO/RAO: -225°/225°, error ±1° C-arm depth:780mm /C-arm opening: 820mm SID:1050mm /Control device Hand switch: Two steps hand switch Foot switch: Dual pedal Console box: 4 joysticks and 5 buttons Image acquisition & processing system PC, CPU2.8GHz, 8GB memory, 1T hard disk, 120GSSD, 2 *USB2.0*1 2 interface touch-screen medical monitors on trolley '19*2 touch-screen monitor on C-arm gantry '13.3 *1 Monitor: DICOM calibration function, meets DICOM standard Fluoroscopy/ جهاز مستعمل ياباني او الماني لا .FLU modes:normal FLU, low-dose FLU, high-dose FLU 3 تقل سنة الصنع عن 2017 -Supporting continuous fluoroscopy and pulse fluoroscopy بالاضافة لبدلات رصاصية واقية من Last image hold function after fluoroscopy C-ARM X-Ray machine الاشعة عدد 3 + واقيات درق Radiology .CINE mode supports 6 different frame: 1,2, 4, 8, 16, 30fps رصاصية عدد 3 - تشمل التكلفة Spot film supports 0.4-320mAs, R' 10 اجور التركيب + طاولة مخصصة Loading time in spot film mode: 40-4000ms, R' 10 للمريض غير نفوذة للاشعة Image post-processing Window width/level adjustment Recursive noise reduction /Imaging zoom /Magnifier Image rotation Horizontal image Vertical image .Annotation, including polygon, circle, rectangle, line, arrow and test .Line segment and angle measurement function Image data management .New case, worklist research, add new examination to exist patients .Deleting cases, delete one or more examinations of patient ,Researching cases according to registration date, patient name & ID including patients who have completed the examination or who have .not completed the examination Support DICOM 3.0 standard worklist inquiry service, download cases .information from HIS/PACS Image browsing Real-time image browsing supports frame rate adjustment of .playback .Real-time image browsing supports current sequence loop play Real-time browsing supports thumbnails to browse all sequences of .current cases Print & storage DICOM film printing is supported the processed DICOM images can .be transferred to a DICOM printer for printing into film Print device management: Multiple printers can be supported at the same time and are compatible with a variety of DICOM format .printers Print format: Supports a variety of film sizes, with a variety of vertical .and horizontal layout functions Support the images export of selected cases to USB flash drives and .CD/DVD discs Support the images export of selected cases to USB flash drives and .CD/DVD discs

CRA/CAU: 90°/45°